



# SCGB2A1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA020816LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O75556
<b>Immunogen</b>	Recombinant Human Mammaglobin-B protein (19-95AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Mammaglobin-B (Lacryglobin) (Lipophilin-C) (Mammaglobin-2) (Secretoglobin family 2A member 1), SCGB2A1, LIPHC MGB2 UGB3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	SCGB2A1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.