



SCGB2A1 Antibody, Biotin conjugated

Product Code	CSB-PA020816LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75556
Immunogen	Recombinant Human Mammaglobin-B protein (19-95AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mammaglobin-B (Lacryglobin) (Lipophilin-C) (Mammaglobin-2) (Secretoglobin family 2A member 1), SCGB2A1, LIPHC MGB2 UGB3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SCGB2A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.