



SCGN Antibody, HRP conjugated

Product Code	CSB-PA020821LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76038
Immunogen	Recombinant Human Secretagogen protein (1-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretagogen, SCGN, SECRET
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCGN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.