



SCGN Antibody, Biotin conjugated

Product Code	CSB-PA020821LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76038
Immunogen	Recombinant Human Secretagogin protein (1-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretagogin, SCGN, SECRET
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCGN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.