



SCLT1 Antibody, HRP conjugated

Product Code	CSB-PA020826LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NL6
Immunogen	Recombinant Human Sodium channel and clathrin linker 1 protein (1-79AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium channel and clathrin linker 1 (Sodium channel-associated protein 1), SCLT1, SAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SCLT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.