



SDSL Antibody, FITC conjugated

Product Code	CSB-PA020927LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GA7
Immunogen	Recombinant Human Serine dehydratase-like protein (39-169AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine dehydratase-like (L-serine deaminase) (L-serine dehydratase/L-threonine deaminase) (L-threonine dehydratase) (TDH) (EC 4.3.1.19) (Serine dehydratase 2) (SDH 2) (EC 4.3.1.17), SDSL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SDSL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.