



SEC22B Antibody, Biotin conjugated

Product Code	CSB-PA020945LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75396
Immunogen	Recombinant Human Vesicle-trafficking protein SEC22b protein (2-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle-trafficking protein SEC22b (ER-Golgi SNARE of 24 kDa) (ERS-24) (ERS24) (SEC22 vesicle-trafficking protein homolog B) (SEC22 vesicle-trafficking protein-like 1), SEC22B, SEC22L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC22B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.