



SELE Antibody, HRP conjugated

Product Code	CSB-PA020975EB01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P98110
Immunogen	Recombinant Pig E-selectin protein (23-429AA)
Raised In	Rabbit
Species Reactivity	Pig
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E-selectin (CD62 antigen-like family member E) (Endothelial leukocyte adhesion molecule 1) (ELAM-1) (Leukocyte-endothelial cell adhesion molecule 2) (LECAM2) (CD antigen CD62E), SELE
Immunogen Species	Pig
Research Area	Others
Target Names	SELE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.