



SLAMF8 Antibody, HRP conjugated

Product Code	CSB-PA021370LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0V8
Immunogen	Recombinant Human SLAM family member 8 protein (23-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLAM family member 8 (B-lymphocyte activator macrophage expressed) (BCM-like membrane protein) (CD antigen CD353), SLAMF8, BLAME
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SLAMF8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.