



SLAMF8 Antibody, HRP conjugated

Product Code	CSB-PA021370LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0V8
Immunogen	Recombinant Human SLAM family member 8 protein (23-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLAM family member 8 (B-lymphocyte activator macrophage expressed) (BCM-like membrane protein) (CD antigen CD353), SLAMF8, BLAME
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SLAMF8