



SLAMF8 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA021370LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9P0V8 |
| Immunogen | Recombinant Human SLAM family member 8 protein (23-233AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | SLAM family member 8 (B-lymphocyte activator macrophage expressed) (BCM-like membrane protein) (CD antigen CD353), SLAMF8, BLAME |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | SLAMF8 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |