



SLC17A3 Antibody, Biotin conjugated

Product Code	CSB-PA021419LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00476
Immunogen	Recombinant Human Sodium-dependent phosphate transport protein 4 protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent phosphate transport protein 4 (Na(+)/PI cotransporter 4) (Sodium/phosphate cotransporter 4) (Solute carrier family 17 member 3), SLC17A3, NPT4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC17A3