



SLC25A14 Antibody, HRP conjugated

Product Code	CSB-PA021480LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95258
Immunogen	Recombinant Human Brain mitochondrial carrier protein 1 protein (55-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Brain mitochondrial carrier protein 1 (BMCP-1) (Mitochondrial uncoupling protein 5) (UCP 5) (Solute carrier family 25 member 14), SLC25A14, BMCP1 UCP5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.