



SLC28A1 Antibody, FITC conjugated

Product Code	CSB-PA021539LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00337
Immunogen	Recombinant Human Sodium/nucleoside cotransporter 1 protein (1-80AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/nucleoside cotransporter 1 (Concentrative nucleoside transporter 1) (CNT 1) (hCNT1) (Na(+)/nucleoside cotransporter 1) (Sodium-coupled nucleoside transporter 1) (Solute carrier family 28 member 1), SLC28A1, CNT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC28A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.