



# SLC31A2 Antibody, HRP conjugated

Product Code	CSB-PA021577LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15432
Immunogen	Recombinant Human Probable low affinity copper uptake protein 2 protein (44-93AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable low affinity copper uptake protein 2 (Copper transporter 2) (hCTR2) (Solute carrier family 31 member 2), SLC31A2, COPT2 CTR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC31A2