



SLC31A2 Antibody, Biotin conjugated

Product Code	CSB-PA021577LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15432
Immunogen	Recombinant Human Probable low affinity copper uptake protein 2 protein (44-93AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable low affinity copper uptake protein 2 (Copper transporter 2) (hCTR2) (Solute carrier family 31 member 2), SLC31A2, COPT2 CTR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC31A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.