



SLC35B1 Antibody, HRP conjugated

Product Code	CSB-PA021588LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78383
Immunogen	Recombinant Human Solute carrier family 35 member B1 protein (106-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 35 member B1 (UDP-galactose transporter-related protein 1) (UGTrel1) (hUGTrel1), SLC35B1, UGTREL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC35B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.