



SLC5A10 Antibody, Biotin conjugated

Product Code	CSB-PA021676LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0PJK1
Immunogen	Recombinant Human Sodium/glucose cotransporter 5 protein (239-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/glucose cotransporter 5 (Na+)/glucose cotransporter 5) (Solute carrier family 5 member 10), SLC5A10, SGLT5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC5A10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.