



SLC5A2 Antibody, HRP conjugated

Product Code	CSB-PA021679LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31639
Immunogen	Recombinant Human Sodium/glucose cotransporter 2 protein (335-422AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/glucose cotransporter 2 (Na(+)/glucose cotransporter 2) (Low affinity sodium-glucose cotransporter) (Solute carrier family 5 member 2), SLC5A2, SGLT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC5A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.