



SLC6A6 Antibody, FITC conjugated

Product Code	CSB-PA021704LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31641
Immunogen	Recombinant Human Sodium- and chloride-dependent taurine transporter protein (1-49AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC6A6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.