



SLITRK3 Antibody, HRP conjugated

Product Code	CSB-PA021772LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94933
Immunogen	Recombinant Human SLIT and NTRK-like protein 3 protein (746-977AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLIT and NTRK-like protein 3, SLITRK3, KIAA0848
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLITRK3