





SLURP1 Antibody

Product Code	CSB-PA021784LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55000
Immunogen	Recombinant Human Secreted Ly-6/uPAR-related protein 1 protein (23-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted Ly-6/uPAR-related protein 1 (SLURP-1) (ARS component B) (ARS(component B)-81/S) (Anti-neoplastic urinary protein) (ANUP), SLURP1, ARS
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SLURP1

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA021784LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA021784LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.