



SLURP1 Antibody, Biotin conjugated

Product Code	CSB-PA021784LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55000
Immunogen	Recombinant Human Secreted Ly-6/uPAR-related protein 1 protein (23-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted Ly-6/uPAR-related protein 1 (SLURP-1) (ARS component B) (ARS(component B)-81/S) (Anti-neoplastic urinary protein) (ANUP), SLURP1, ARS
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SLURP1