





SOCS3 Antibody

Product Code	CSB-PA022392EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14543
Immunogen	Recombinant Human Suppressor of cytokine signaling 3 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppressor of cytokine signaling 3 (SOCS-3) (Cytokine-inducible SH2 protein 3) (CIS-3) (STAT-induced STAT inhibitor 3) (SSI-3), SOCS3, CIS3 SSI3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SOCS3

Image



Immunocytochemistry analysis of human liver cancer using CSB-PA022392EA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA022392EA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.