



SOCS3 Antibody, HRP conjugated

Product Code	CSB-PA022392EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14543
Immunogen	Recombinant Human Suppressor of cytokine signaling 3 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Suppressor of cytokine signaling 3 (SOCS-3) (Cytokine-inducible SH2 protein 3) (CIS-3) (STAT-induced STAT inhibitor 3) (SSI-3), SOCS3, CIS3 SSI3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	SOCS3