




SPRR1B Antibody

Product Code	CSB-PA022611LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22528
Immunogen	Recombinant Human Cornifin-B protein (1-89AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cornifin-B (14.9 kDa pancornulin) (Small proline-rich protein IB) (SPR-IB), SPRR1B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPRR1B

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA022611LA01HU at dilution of 1:100