



SPRR1B Antibody, FITC conjugated

Product Code	CSB-PA022611LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22528
Immunogen	Recombinant Human Cornifin-B protein (1-89AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cornifin-B (14.9 kDa pancornulin) (Small proline-rich protein IB) (SPR-IB), SPRR1B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPRR1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.