



SRP54 Antibody, FITC conjugated

Product Code	CSB-PA022675LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61011
Immunogen	Recombinant Human Signal recognition particle 54 kDa protein (1-504AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal recognition particle 54 kDa protein (SRP54), SRP54
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SRP54
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.