



SRP54 Antibody, FITC conjugated

Product Code	CSB-PA022675LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61011
Immunogen	Recombinant Human Signal recognition particle 54 kDa protein (1-504AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal recognition particle 54 kDa protein (SRP54), SRP54
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SRP54