





SST Antibody

Product Code	CSB-PA022723OJ01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61278
Immunogen	Peptide sequence from Human Somatostatin protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Somatostatin (Growth hormone release-inhibiting factor) [Cleaved into: Somatostatin-28; Somatostatin-14], SST
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SST

Image



Immunohistochemistry analysis of human pancreatic cancer using CSB-PA022723OJ01HU at dilution of 1:100



Immunohistochemistry analysis of human colon cancer using CSB-PA022723OJ01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.