



ST8SIA5 Antibody, HRP conjugated

Product Code	CSB-PA022778LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15466
Immunogen	Recombinant Human Alpha-2,8-sialyltransferase 8E protein (101-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2,8-sialyltransferase 8E (EC 2.4.99.-) (Sialyltransferase 8E) (SIAT8-E) (Sialyltransferase St8Sia V) (ST8SiaV), ST8SIA5, SIAT8E
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ST8SIA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.