



ST8SIA5 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA022778LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O15466 |
| Immunogen | Recombinant Human Alpha-2,8-sialyltransferase 8E protein (101-376AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Alpha-2,8-sialyltransferase 8E (EC 2.4.99.-) (Sialyltransferase 8E) (SIAT8-E) (Sialyltransferase St8Sia V) (ST8SiaV), ST8SIA5, SIAT8E |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ST8SIA5 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |