



STAT1 Antibody, Biotin conjugated

Product Code	CSB-PA022810LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42224
Immunogen	Recombinant Human Signal transducer and activator of transcription 1-alpha/beta protein (63-196AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal transducer and activator of transcription 1-alpha/beta (Transcription factor ISGF-3 components p91/p84), STAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.