



STK17B Antibody, FITC conjugated

Product Code	CSB-PA022838LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94768
Immunogen	Recombinant Human Serine/threonine-protein kinase 17B protein (235-366AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 17B (EC 2.7.11.1) (DAP kinase-related apoptosis-inducing protein kinase 2), STK17B, DRAK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	STK17B