



STK17B Antibody, Biotin conjugated

Product Code	CSB-PA022838LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94768
Immunogen	Recombinant Human Serine/threonine-protein kinase 17B protein (235-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase 17B (EC 2.7.11.1) (DAP kinase-related apoptosis-inducing protein kinase 2), STK17B, DRAK2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	STK17B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.