




STT3A Antibody

Product Code	CSB-PA022885LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46977
Immunogen	Recombinant Human Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3A protein (39-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3A (Oligosaccharyl transferase subunit STT3A) (STT3-A) (EC 2.4.99.18) (B5) (Integral membrane protein 1) (Transmembrane protein TMC), STT3A, ITM1 TMC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	STT3A

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA022885LA01HU at dilution of 1:100