



STT3A Antibody, Biotin conjugated

Product Code	CSB-PA022885LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46977
Immunogen	Recombinant Human Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3A protein (39-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3A (Oligosaccharyl transferase subunit STT3A) (STT3-A) (EC 2.4.99.18) (B5) (Integral membrane protein 1) (Transmembrane protein TMC), STT3A, ITM1 TMC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	STT3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.