



SYCN Antibody, HRP conjugated

Product Code	CSB-PA022989LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0VAF6
Immunogen	Recombinant Human Syncollin protein (22-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syncollin (Insulin synthesis-associated protein 1), SYCN, INSSA1 SYL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SYCN