





TAS2R41 Antibody

Product Code	CSB-PA023159LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59536
Immunogen	Recombinant Human Taste receptor type 2 member 41 protein (205-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 41 (T2R41) (Taste receptor type 2 member 59) (T2R59), TAS2R41
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TAS2R41

Image



Immunofluorescent analysis of MCF-7 cells using CSB-PA023159LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA023159LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.