



TAS2R41 Antibody, Biotin conjugated

Product Code	CSB-PA023159LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59536
Immunogen	Recombinant Human Taste receptor type 2 member 41 protein (205-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 41 (T2R41) (Taste receptor type 2 member 59) (T2R59), TAS2R41
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TAS2R41
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.