



TAZ Antibody, FITC conjugated

Product Code	CSB-PA023181EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F6Y2X3
Immunogen	Recombinant Human Tafazzin family protein (95-262AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tafazzin family protein, TAZ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.