



Tbx15 Antibody

Product Code	CSB-PA023247LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O70306
Immunogen	Recombinant Mouse T-box transcription factor TBX15 protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-box transcription factor TBX15 (T-box protein 15) (MmTBx8) (T-box transcription factor TBX14) (T-box protein 14), Tbx15, Tbx14 Tbx8
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Tbx15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.