



Tbx15 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA023247LD01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O70306 |
| Immunogen | Recombinant Mouse T-box transcription factor TBX15 protein (1-120AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | T-box transcription factor TBX15 (T-box protein 15) (MmTBx8) (T-box transcription factor TBX14) (T-box protein 14), Tbx15, Tbx14 Tbx8 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Others |
| Target Names | Tbx15 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |