



TCN2 Antibody, HRP conjugated

Product Code	CSB-PA023318LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20062
Immunogen	Recombinant Human Transcobalamin-2 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcobalamin-2 (TC-2) (Transcobalamin II) (TC II) (TCII), TCN2, TC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TCN2