



TCTA Antibody, Biotin conjugated

Product Code	CSB-PA023331LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57738
Immunogen	Recombinant Human T-cell leukemia translocation-altered gene protein (62-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell leukemia translocation-altered gene protein (T-cell leukemia translocation-associated gene protein), TCTA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TCTA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.