




tgfb1 Antibody

Product Code	CSB-PA023446LA01OEI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O93449
Immunogen	Recombinant Oncorhynchus mykiss Transforming growth factor beta-1 proprotein protein (271-382AA)
Raised In	Rabbit
Species Reactivity	Oncorhynchus mykiss, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1 proprotein [Cleaved into: Latency-associated peptide (LAP)], (Transforming growth factor beta-1) (TGF-beta-1) tgfb1
Immunogen Species	Oncorhynchus mykiss
Research Area	Others
Target Names	tgfb1

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, A549 whole cell lysate, HepG2 whole cell lysate

All lanes: tgfb1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.