



tgfb1 Antibody, HRP conjugated

Product Code	CSB-PA023446LB01OEI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O93449
Immunogen	Recombinant Oncorhynchus mykiss Transforming growth factor beta-1 proprotein protein (271-382AA)
Raised In	Rabbit
Species Reactivity	Oncorhynchus mykiss
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1 proprotein [Cleaved into: Latency-associated peptide (LAP)], (Transforming growth factor beta-1) (TGF-beta-1) tgfb1
Immunogen Species	Oncorhynchus mykiss
Research Area	Others
Target Names	tgfb1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.