



# Timp1 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA023560LD01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P12032
<b>Immunogen</b>	Recombinant Mouse Metalloproteinase inhibitor 1 protein (25-205AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Metalloproteinase inhibitor 1 (Collagenase inhibitor 16C8 fibroblast) (Erythroid-potentiating activity) (EPA) (TPA-S1) (TPA-induced protein) (Tissue inhibitor of metalloproteinases 1) (TIMP-1), Timp1, Timp Timp-1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	Timp1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.