



# TK Antibody, HRP conjugated

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA023577HB01HWY                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | P03176                                                                           |
| Immunogen           | Recombinant Human herpesvirus 1 Thymidine kinase protein (1-376AA)               |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Human herpesvirus 1                                                              |
| Tested Applications | ELISA                                                                            |
| Form                | Liquid                                                                           |
| Conjugate           | HRP                                                                              |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified                                                         |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Alias               | Thymidine kinase (EC 2.7.1.21), TK                                               |
| Immunogen Species   | Human herpesvirus 1                                                              |
| Research Area       | Others                                                                           |
| Target Names        | TK                                                                               |