



TK Antibody, HRP conjugated

Product Code	CSB-PA023577HB01HWY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03176
Immunogen	Recombinant Human herpesvirus 1 Thymidine kinase protein (1-376AA)
Raised In	Rabbit
Species Reactivity	Human herpesvirus 1
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymidine kinase (EC 2.7.1.21), TK
Immunogen Species	Human herpesvirus 1
Research Area	Others
Target Names	TK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.